


The Features of the Easy-PREPCCC instrument

KUTUWA Sangyo co.,LTD.

Features

- ① Side opened during maintenance)
 - ② Cartridge type transfer tube
 - ③ Flange type cartridge column
 - ④ With hanging rotor motor
 - ⑤ Touch panel
- Stop time setting operation (PUMP stop signal)
 - Solvent front alarm
 - Peak detection valve switching contact output
 - Leak sensor (option)
 - Alarm stop pump, motor stop

① Side door (opened during maintenance)




② Cartridge type transfer tube can be attached and removed easily.




② Cartridge type transfer tube can be attached and removed easily.



Cartridge type transfer tube can be fixed with screws.




③ Flange type cartridge column can be connected easily with liquid leakage prevention.




Flange type screw

Flange type cartridge column.

④ With hanging rotor motor



⑤ Touch panel



- Stop time setting operation (PUMP stop signal)
- Solvent front alarm
- Peak detection valve switching contact output
- Leak sensor (option)
- Alarm stop pump, motor stop